

The slide features three decorative curved lines in the corners, each composed of multiple overlapping layers in shades of light green and blue. One arc is in the top right, another in the bottom left, and a third in the bottom right.


# Assignment

Knee Preop Planning

# Step-1











- Load STL in the web application scene
  - Femur.stl
  - Tibia.stl

*Note: The position of object with respect to world coordinate should not change while importing in the web application.*



## Step-2


- Create Landmarks(Points)


- Femur Center 
- Hip Center 
- Femur Proximal Canal 
- Femur Distal Canal 
- Medial Epicondyle 
- Lateral Epicondyle 
- Distal Medial Pt 
- Distal Lateral Pt 
- Posterior Medial Pt 
- Posterior Lateral Pt 

- Default all buttons are in inactive mode,
- First user clicks on Femur Center Landmark Active Button(so it will change color from light gray to black) then Landmark(Point) creation functionality will activate
- Then the user click on femur bone surface then landmark will be created, also translation (in X-Y-Z) Control will be activated
- If the user want, then he can move the landmark
- Then user either click on Another landmark or Femur Center Landmark Active Button
- So Translation control will not appear in scene, also button will change color to gray
- Then user click on Hip Center Landmark Active button and again Landmark creation functionality will activate and user click on surface and landmark will be created, and translation control will be applied
- Similar way user will create all the 10 landmarks on Femur Bone
- Landmarks position is shown in next page










 Active Mode



 Inactive Mode

## Step-2 (As Result)







## Step-3

- Update Button:

- After all landmarks creation user will click on “Update Button”
- So the first Axis (lines) will create.


- Line Create – Between Femur Center & Hip Center Landmark (Mechanical Axis)
- Line Create – Between Femur Proximal Canal & Femur Distal Canal (Anatomical Axis)
- Line Create – Between Medial Epicondyle & Lateral Epicondyle (TEA-Trans epicondyle Axis)
- Line Create – Between Posterior Medial Pt & Posterior Lateral Pt (PCA- Posterior Condyle Axis)



## Step-4 (Continue)




- Plane Create – Perpendicular to Mechanical Axis




- Project "TEA" Axis on Perpendicular Plane

## Step-4 (Continue)




- Create a line on a perpendicular plane, that line should start from Femur Center Landmark & have 10mm distance Anterior Side , Also this line should be perpendicular to the Projected TEA Axis




- Duplicate Perpendicular plane and make functionality like
- If user click on Positive or Negative button, Then New Plane should be Rotate along to Anterior Line (Anterior Line is a rotation axis for new plane) we can call new plane as varus/valgus plane



## Step-4 (Continue)




- Project Anterior Line on Varus/Valgus Plane
- Then create a on Varus/Valgus Plane, it should be perpendicular to the Projected Anterior line and start from femur center & end at 10mm Lateral side




- Duplicate Varus/Valgus plane and make functionality like
- If user click on Positive or Negative button Then New Plane(Flexion/Extension Plane) should be Rotate along to Lateral Line (Lateral Line is a rotation axis for a new plane) we can call new plane as Flexion/Extension plane

## Step-4 (Continue)




- Again create new plane
- It is parallel plane to flexion/extension plane
- Also it is passing through Distal Medial Pt
- New plane called as distal medial plane



- Again new plane create
- It is parallel to distal medial plane & have 10mm distance
- 10mm distance should be on proximal direction.
- New plane called as distal resection plane


## Step-4 (Continue)



- Resect distal part of femur with help of Distal Resection Plane
- Also, resection can be controlled by a toggle button
  - If toggle is on = Resection will be shown in a scene
  - If toggle is off = Resection will not show in a scene

## Step-4 (Continue)

- Show measurement in a scene
  - Distal Medial = Distance between Distal resection plane and Distal Medial Pt
  - Distal Lateral = Distance between Distal resection plane and Distal Lateral Pt



# Application Feature

- 1 Scene
- Femur and Tibia Bone visible
- 4 Axis & 10 landmarks will be visible
- Plane should be hide
- Landmark creation and modification buttons (10 buttons)
- Update button
- Varus/Valgus Button = Increment and Decrement (°)
- Flexion/Extension Button = Increment and decrement (°)
- Distal Medial Resection= Increment and decrement (mm)
- Toggle button = to control resection visibility

Note: If user modify anything then all child axis and plane will be update  
as example ,  
If Varus/valgus of 3° is changed to 1° then,  
Flexion/extension plane will change, then  
distal medial & distal resection also change  
so resection of femur also change ,  
measurement also change.



A decorative wavy line in the top left corner, transitioning from light green to light blue.

Thank You

A large decorative circular arc on the right side, transitioning from light blue to light green.